

TITANIUM BASE

For implant supported individual hybrid abutments



BENEFITS:

- engaged
- for individual abutment design
- large glue surface, so durable solid connection to the base
- conventional usable and also for CAD/CAM process



APPLICATION:

- for individual glued ceramic abutments (zirconia, lithiumdisilicate) or synthetics (PMMA, PEEK, composite)



TECHNICAL DETAILS:

- Titanium Grade 5 ELI



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures. Products indicated with ★ are available in the USA 

TABLE OF CONTENTS

NT CATALOGUE 2017-1

1	TITANIUM BASE ABUTMENT	10 – 49
2	2-CONnect® ABUTMENT	50 – 117
3	NT-PREFORMS® ABUTMENT	118 – 147
4	DIM ANALOG	148 – 161
5	SURGICAL GUIDE SLEEVES	163 – 165
6	TOOLS & DISCS	167 – 181
7	NT-iDent	182 – 189
8	SPECIFICATION	190 – 191

BEG-SERIES – NEW

is compatible to: BEGO Implant Systems® / Semados®
S, RS-RX, RI-Line - - - - - **12**

BIO-SERIES – NEW

is compatible to: BioHorizons® / internal- - - - - **13**

BRE-SERIES – NEW

is compatible to: Bredent Medical® / Sky® - - - - - **14**

CAM-SERIES

is compatible to: Altatec® / Camlog® - - - - - **15**

COL-SERIES

is compatible to: Altatec® / Conelog® - - - - - **16**

E-SERIES

is compatible to: Nobel Biocare™ / Replace Select®- - - - - **17**

F-SERIES

is compatible to: Nobel Biocare™ / Nobel Active™ - - - - - **21**

H-SERIES

is compatible to: Biomet 3i® / Certain® - - - - - **22**

HIO-SERIES – NEW

is compatible to: HIOSSSEN® / OSSTEM TS® - - - - - **23**

I-SERIES

is compatible to: Biomet 3i® / Osseotite® External Hex® - - **24**

K-SERIES

is compatible to: Nobel Biocare™ / Brånemark® - - - - - **25**

L-SERIES

is compatible to: Straumann / Bone Level®- - - - - **26**

MEX-SERIES – NEW

is compatible to: Medentis Medical® / ICX®- - - - - **29**

MIC-SERIES – NEW

is compatible to: MIS® / C1® - - - - - **30**

MIS-SERIES – NEW

is compatible to: MIS® / SEVEN®- - - - - **31**

N-SERIES

is compatible to: Straumann / SynOcta® - - - - - **32**

NEO-SERIES – NEW

is compatible to: NEOSS Implant System® / NEOSS™ - - - - - **35**

R-SERIES

is compatible to: Zimmer Dental / Tapered Screw-Vent® - - **36**

S-SERIES

is compatible to: DENTSPLY Implants® /
Astra Tech Implant System OsseoSpeed® TX - - - - - **39**

SEV-SERIES – NEW

is compatible to: DENTSPLY Implants® /
Astra Tech Implant System™ EV - - - - - **42**

T-SERIES

is compatible to: DENTSPLY Implants® / Xive®- - - - - **43**

THO-SERIES – NEW

is compatible to: Thommen Medical® / SPI® - - - - - **46**

UPR-SERIES – NEW

is compatible to: Sweden & Martina Implantology™ /
Premium Kohno®- - - - - **47**

UTG-SERIES – NEW

is compatible to: Sweden & Martina Implantology™ /
Premium TG®- - - - - **48**



BEG-SERIES

compatible to: BEGO Implant Systems® / Semados® S, RS-RX, RI-Line



PLATFORM

3,25 / 3,75 mm

4,1 mm

4,5 mm

5,5 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



BEG 80323700

BEG 80410000

BEG 80450000

BEG 80550000

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



BEG 9.S3D3.237

BEG 9.S3D4.100

BEG 9.S3D4.500

BEG 9.S3D5.500

★

NT-MILL LAB ANALOG

Stainless Steel



BEG 5.La3.237

BEG 5.La4.100

BEG 5.La4.500

BEG 5.La5.500

★

IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



BEG TR-00013.2

BEG TR-00024.1

BEG TR-00034.5

BEG TR-00045.5

ABUTMENT SCREW

Recommended tightening
torques 30 Ncm
Titanium Grade 5 ELI



BEG 603255

BEG 603255

BEG 603255

BEG 603255

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures. Products indicated with ★ are available in the USA

BIO-SERIES

compatible to: BioHorizons® / internal



PLATFORM

3,0 mm

3,5 mm

4,5 mm

5,7 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



BIO 8030000

BIO 8035000

BIO 8045000

BIO 8057000

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



BIO 9.S3D3.000

BIO 9.S3D3.500

BIO 9.S3D4.500

BIO 9.S3D5.700

★

NT-MILL LAB ANALOG

Stainless Steel



BIO 5.La3.000

BIO 5.La3.500

BIO 5.La4.500

BIO 5.La5.700

★

IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



BIO TR-00013.0

BIO TR-00023.5

BIO TR-00034.5

BIO TR-00045.7

ABUTMENT SCREW

Recommended tightening
torques 30 Ncm
Titanium Grade 5 ELI



BIO 603057

BIO 603057

BIO 603057

BIO 603057

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

BRE-SERIES

compatible to: Bredent Medical® / Sky®



PLATFORM

3,5 mm

4,0 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



BRE 80350000

BRE 80400000

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



BRE 9.S3D3.500

BRE 9.S3D4.000

★

NT-MILL LAB ANALOG

Stainless Steel



BRE 5.La3.500

BRE 5.La4.000

★

IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



BRE TR-00013.5

BRE TR-00024.0

ABUTMENT SCREW

Recommended tightening
torques 25 Ncm
Titanium Grade 5 ELI



BRE 603540

BRE 603540

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

CAM-SERIES

compatible to: Altatec® / Camlog®



PLATFORM

3,3 mm

3,8 mm

4,3 mm

5,0 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screws
Titanium Grade 5 ELI



CAM 80330000

CAM 80380000

CAM 80430000

CAM 80500000



SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



CAM 9.S3D3.300

CAM 9.S3D3.800

CAM 9.S3D4.300

CAM 9.S3D5.060



NT-MILL LAB ANALOG

Stainless Steel



CAM 5.La3.300

CAM 5.La3.800

CAM 5.La4.300

CAM 5.La5.060



IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



CAM TR-00013.3

CAM TR-00023.8

CAM TR-00034.3

—

LAB SCREW

coloured screw



CAM 6233843

CAM 6233843

CAM 6233843

CAM 6250600



ABUTMENT SCREW

Recommended tightening
torques 20 Ncm
Titanium Grade 5 ELI



CAM 6033843

CAM 6033843

CAM 6033843

CAM 6150600



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

COL-SERIES

compatible to: Altatec® / Conelog®



PLATFORM

3,3 mm

3,8 / 4,3 mm

5,0 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screws
Titanium Grade 5 ELI



Gingiva Height 0,8 mm
..... 2,0 mm

COL 80330008

COL 80384308

COL 80500008

COL 80330020

COL 80384320

COL 80500020

SCAN BODY

for Titanium Base
incl. coloured lab screw
PEEK



COL 703300

COL 703843

COL 705000

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



COL 9.S3D3.300

COL 9.S3D3.843

COL 9.S3D5.000

★

NT-MILL LAB ANALOG

Stainless Steel



COL 5.La3.300

COL 5.La3.843

COL 5.La5.000

★

IMPLANT PICK-UP

for open-tray with
integrated screw
Coming soon



COL TR-00013.3

COL TR-00023.8

—

LAB SCREW

coloured screw



COL 623343

COL 623343

COL 625000

ABUTMENT SCREW

Recommended tightening
torques 20 Ncm
Titanium Grade 5 ELI



COL 603343

COL 603343

COL 605000

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

E-SERIES

compatible to: Nobel Biocare™ / Replace Select®



PLATFORM

3,5 mm NP

4,3 mm RP

5,0 mm WP

6,0 mm

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



E 800

E 810

E 820

E 830

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



E 9.S3D3.500

E 9.S3D4.300

E 9.S3D5.000

E 9.S3D6.000

★

LAB ANALOG

Stainless Steel



E 50

E 51

E 52

E 53

★

IMPLANT PICK-UP

for open-tray with
integrated screw



E TR-NP013.5

E TR-RP024.3

E TR-WP035.0

—

ABUTMENT SCREW

Recommended tightening
torques 35 Ncm
Titanium Grade 5 ELI



E 60

E 61

E 61

E 61

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 

E-SERIES

compatible to: Nobel Biocare™ / Replace Select®



PLATFORM

3,5 mm NP

4,3 mm RP

5,0 mm WP

STRAIGHT ABUTMENT

Gingiva Height 1,0 mm
incl. screw
Titanium Grade 5 ELI



E 100-1

E 110-1

E 120-1

★

STRAIGHT ABUTMENT

Gingiva Height 2,5 mm
incl. screw
Titanium Grade 5 ELI



E 100

E 110

E 120

★

UCLA ABUTMENT

engaging
incl. screw



E 11.CA3.500

E 11.CA4.300

E 11.CA5.000

ABUTMENT SCREW

Recommended tightening
torques 35 Ncm
Titanium Grade 5 ELI



E 60

E 61

E 61

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

E-SERIES

compatible to: Nobel Biocare™ / Replace Select®



PLATFORM

3,5 mm NP

4,3 mm RP

5,0 mm WP

ANGLED ABUTMENT 16°

angled over surface
Gingiva Height 1,0 mm
incl. screw
Titanium Grade 5 ELI



E 200-1-1

E 210-1-1

E 220-1-1

ANGLED ABUTMENT 16°

angled over edge
Gingiva Height 1,0 mm
incl. screw
Titanium Grade 5 ELI



E 200-2-1

E 210-2-1

E 220-2-1

ANGLED ABUTMENT 16°

angled over surface
Gingiva Height 2,5 mm
incl. screw
Titanium Grade 5 ELI



E 200-1

E 210-1

E 220-1

ANGLED ABUTMENT 16°

angled over edge
Gingiva Height 2,5 mm
incl. screw
Titanium Grade 5 ELI



E 200-2

E 210-2

E 220-2

ABUTMENT SCREW

Recommended tightening
torques 35 Ncm
Titanium Grade 5 ELI



E 60

E 61

E 61

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA